

Tapestri Single-Cell DNA and DNA + Protein Sequencing User Guide [Workflow Changes](#)

	Chapter	Step	V2	V3
0	Firmware	Instrument Firmware		<a href="#">Update/ensure V3 compatible firmware</a>
1	Prepare Cell Suspension	Dilute Cell Suspension	<b>3,000 - 4,000</b> cells/ $\mu$ L	<b>2,800 - 3,200</b> cells/ $\mu$ L
2	Encapsulate Cells	Reagent Retrieval		Equilibrate Encapsulation Oil to RT <b>from 4°C storage</b>
		Prepare Lysis Mix	Prepare 100 $\mu$ L Lysis Mix, <b>use 90 <math>\mu</math>L</b>	Prepare 70.1 $\mu$ L Lysis Mix, <b>use 60 <math>\mu</math>L</b>
		Post Encapsulation	<b>Remove most oil</b> leaving <5 $\mu$ L only	Remove oil until <b>total volume (emulsion + oil) = 100 <math>\mu</math>L</b>
3	Lyse & Digest Cells	No Change		
4	Barcode Cells	Reagent Retrieval	Equilibrate Barcoding Beads to RT <b>from 4°C storage</b>	Equilibrate Barcoding Beads to RT <b>from -20°C storage</b> Equilibrate Barcoding Oil to RT <b>from 4°C storage</b>
		Barcoding Bead (BC) Preparation	<b>Use Barcoding Beads as is</b>	<b>Add 67 <math>\mu</math>L of Barcode Mix to Barcoding Beads</b> prior to vortexing
		Barcoding	Load <b>200 <math>\mu</math>L Barcoding Beads</b> into reservoir 7	Load <b>250 <math>\mu</math>L Barcoding Beads</b> into reservoir 7
			Load <b>250 <math>\mu</math>L of Barcode Mix</b> into reservoir 8	Load <b>200 <math>\mu</math>L of Barcode Mix</b> into reservoir 8
Post Barcoding	Remove <b>up to 120 <math>\mu</math>L oil</b> from each tube	Remove oil until the <b>total volume (emulsion + oil) = 100 <math>\mu</math>L</b> from each tube		
5	Targeted Amplification	No Change		
6	Cleanup PCR Products	Digest PCR Product	<b>Pool tubes 1-4 and 5-8 into separate tubes</b> and add 20 $\mu$ L DNA Clean up Buffer and 12 $\mu$ L Clean up Enzyme Incubate in <b>four 0.2 mL tubes</b>	<b>Pool tubes 1-8 into one tube</b> and add 40 $\mu$ L DNA Clean up Buffer and 24 $\mu$ L Clean up Enzyme Incubate in <b>one 1.5 mL tube</b>
		Clean up PCR product		After enzymatic digest, spin down tube and <b>transfer supernatant leaving behind any visible pellet</b>
		Clean up PCR product	Perform AMPure XP bead cleanup in <b>2 tubes</b>	Perform AMPure XP bead cleanup in a <b>single tube</b>
		DNA Library Cleanup I	Perform one 0.72x (DNA) or 0.70x (DNA+Protein) AMPure XP cleanup	Perform one 0.72x (DNA) or 0.70x (DNA+Protein) AMPure XP cleanup <b>followed by a second 0.76x AMPure XP cleanup</b>
7	Library PCR	DNA Library Cleanup II	Perform <b>one 0.69x AMPure XP</b> cleanup	Perform one 0.69x AMPure XP cleanup <b>followed by a second 0.72x AMPure XP cleanup</b>
8	Quantify and Normalize Library	No Change		
9	Sequence Library	Required reads	<b>~5,000 cells expected</b>	<b>~11,000 cells expected</b>