

## Tapestri® Genotype Multiplexing

### IMPORTANT

This protocol describes the necessary steps to multiplex up to three **unrelated** samples upstream of Encapsulation on the Tapestri® instrument. Germline variants are used for sample demultiplexing using the Tapestri Pipeline. For accurate demultiplexing, each sample must comprise at least 0.1% of the pooled cell suspension.

### DNA-only Applications

**NOTE** Follow the Tapestri® Single-Cell DNA v3 User Guide, PN MB05-0017, [Steps 1.1–1.8](#) (through cell resuspension in DPBS, before centrifugation). Perform [Steps 1.2–1.5](#) for each sample one at a time and store on ice, then process all samples in parallel.

1. **Quantify the cells and assess viability** using an automated cell counter or hemocytometer following best practices and the manufacturer's instructions.

**NOTE** The recommended input for cell encapsulation is 100,000 cells (30,000 minimum). It is recommended to use 150,000 cells in [Step 2](#) below (50,000 minimum), to account for cell loss.

2. **Pool cells** from up to three samples for a **total of 150,000 cells** in a new 1.5 mL DNA LoBind tube. If needed, use a 15 mL conical tube to accommodate larger volumes.

**NOTE** If using a diploid cell line as a control for CNV measurement, add enough diploid cells to constitute 5% of the final pool.

3. Centrifuge the multiplexed cell suspension at **400 x g for 5 minutes** at 4 °C in a swinging bucket.
4. **Place cells on ice** and proceed **immediately** to [Step 1.9](#) of the Tapestri® Single-Cell DNA v3 User Guide, PN MB05-0017.

### DNA + Protein Applications

**NOTE** Follow the Tapestri® Single-Cell DNA + Protein v3 User Guide, PN MB05-0018, [Steps 1.1–1.16](#) (through cell resuspension in Cell Staining Buffer). Perform [Steps 1.3–1.9](#) for each sample one at a time and store on ice, then process all samples in parallel.

1. **Quantify the cells and assess viability** using an automated cell counter or hemocytometer following best practices and the manufacturer's instructions.

**NOTE** The recommended input for cell staining is 1M cells (400,000 minimum). It is recommended to use 1.5M cells in [Step 2](#) below (600,000 minimum), to account for cell loss.

2. **Pool cells** from up to three samples for a **total of 1.5M cells** in a new 1.5 mL DNA LoBind tube. If needed, use a 15 mL conical tube to accommodate larger volumes.

**NOTE** If using a diploid cell line as a control for CNV measurement, add enough diploid cells to constitute 5% of the final pool.

3. **Place cells on ice** and proceed **immediately** to [Step 1.17](#) of the Tapestri® Single-Cell DNA + Protein v3 User Guide, PN MB05-0018.