

Tapestri Single-Cell DNA to DNA + RNA Workflow Changes

	Chapter	Step	DNA	DNA + RNA
	Reagents			Do not use Lysis Mix or Cell Buffer from DNA kit
	Firmware	Instrument Firmware	No Change	No Change
	Panel Design	Experimental Design	Panel design in Designer	Panel design supported by Mission Bio
1	Prepare Cell Suspension		Dilute in CB to 3,000 cells/uL	Dilute cells in DPBS to 4,500 cells/uL
2	Encapsulate Cells	Lysis Mix Preparation	Standard Lysis Mix with DNA Rev Primer	RNA Lysis Mix with RNA Enzyme/reagents/RNA Rev Primer
3	Lyse & Digest Cells		1 hr lysis time	2 hrs lysis time
4	Barcode Cells	Barcode Mix Preparation	Barcoding mix with DNA Fwd Primers	Barcoding mix with RNA Fwd and DNA Fwd + Rev primers
5	Targeted Amplification		No Change	No Change (use DNA panel size annealing time)
6	Cleanup PCR Products	DNA Library Cleanup I	0.72X, 0.76X AmpureXP cleanups	0.72X, AmpureXP cleanup, agarose gel cleanup , 0.92X AmpureXP cleanup
7	Library PCR	DNA Library Cleanup II	0.69X, 0.72X AmpureXP cleanups	Side by side 2x 0.69X AmpureXP cleanups of DNA and RNA libraries
8	Quantify and Normalize Library		DNA library with ~460 bp peak	DNA and RNA libraries with ~460 bp peak
9	Sequence Library	Required reads	~80x coverage per cell/amplicon	DNA ~80x coverage per cell/amplicon, RNA 100x per cell/gene